



Double Beam Spectrophotometer
U-2900/2910

HITACHI
Inspire the Next



U-2900/2910

Spectrophotometer simpler to use and higher in dependability thanks to mounting of a large-size color LCD.

High resolution satisfying European Pharmacopoeia (spectral bandpass: 1.5 nm).

Trace-amount measurement in biotechnological field, etc. can be carried out using optional 50, 25 and 5 mL micro-volume cells.

- An abundance of optional accessories for various applications have been lined up.
- Function for automatic data storage in a USB memory is incorporated.
- Validation function and self-diagnostic function are standard with this instrument.
- The UV Solutions control program available at option enables the user to control operations from a personal computer.

(The Model U-2910 is a special model for PC control.)

*: Available at option



Model 124



Model 100-50



Model 150-20



Model U-2000



Model U-2001




Model U-2800



**U-2900
/2910**

Stable optics due to double beam



In double beam design, the energy of the light source is divided into two with a half mirror so that one passes through the reference side, and the other through the sample side, which is unavailable with the single beam design.


Since the reference-side energy is also incident on a detector, photometry is carried out on the basis of this signal. Therefore, an energy change in the light source can be compensated to ensure stable measurement for a long time.

Incorporation of stigmatic concave diffraction grating

The optics of this instrument adopts the Seya-Namioka monochromator widespread as a representative concave diffraction grating monochromator.

Because a concave diffraction grating has both beam condensing and dispersing functions, an optical system can be configured with fewer mirrors. In a spectrophotometer, use of fewer mirrors signifies a shorter optical path, thus giving rise to an aberration-free bright optics.

For elimination of the aberrations which were essentially unavoidable in the past, a stigmatic concave diffraction grating has been developed by applying Hitachi's original technology. As a result, a higher resolution has been realized.




A higher resolution has been achieved by eliminating coma from the Seya-Namioka monochromator which is the most popular concave diffraction grating monochromator. Its grooving is supported by the only ruling engine in Japan. The diffraction gratings of Model U-2900/2910 have also been made with this machine.

Spectral bandpass 1.5 nm satisfying European Pharmacopoeia

European Pharmacopoeia requires a ratio of 1.5 or larger when measuring 0.02% (V/V) solutions of toluene in hexane. In the spectrum at right, a satisfactory value of 1.6 can be confirmed.

(Sample: 0.02% solution of toluene in n-hexane)



Mounting of 26.4 cm color LCD, the largest in this product class



* The size in the photo is 1/2 of the actual LCD size.

Due to a 640×480 color LCD with backlight, working curves, spectra and other minute displays are easier to view, and characters are bigger than before, which ensures a high visibility. A dedicated keyboard (for U-2900 alone) is quite helpful for baseline measurement and other operations.




* The photo size is 1/2 of the actual keyboard size.

Introduction of Measurement Modes

Photometry

The concentration of an unknown sample can be determined in comparison with the known concentration of a standard sample.


The linear/quadratic regression curve or polygonal line approximation curve is drawn by measuring plural standard samples (max. 20 samples) to calculate the concentration. In addition, judgment of the upper and lower concentration limits can be set and quantitative analysis can also be carried out via factor input.



Wavelength scan



A spectrum resulting from scanning over a desired wavelength range within 190 to 1,100 nm can be displayed. Because each substance has a unique spectrum, its characteristics can be examined.

After measurement, peak search, smoothing and other data processing can be effected. And, a repeat scan capability helps the user trace a chemical reaction process. Baseline correction is also available.




Time scan

A spectrum can be drawn by following a change of photo-metric value at a single wavelength with the lapse of time. This function is used for analysis of enzymatic reactions. Enzyme activity is measured according to a change of absorbance within the set time period. Peak detection, smoothing, kinetic assay and other data processing can be performed.



Multiple-wavelength measurement

Purity of a nucleic acid can be calculated using a ratio of absorbance values at 2 wavelengths (A_{260}/A_{280}). Furthermore, measurement can be performed while automatically shifting measuring wavelength to a maximum of 6 wavelengths. This is convenient when only measured data at each wavelength is desired.



Highly dependable system incorporating automatic calibration and self-diagnostic function.

A highly dependable system was designed by incorporating lots of functions for memory check, wavelength drive check, lamp ignition check, automatic wavelength calibration, lamp ignition time display, etc.


(This function is standard with the Model U-2900 main unit. The Model U-2910 provides the equivalent functions via PC control.)

Validation function


This is a function for checking the performance of Model U-2900/2910 main unit.

The following items can be checked.

- Wavelength accuracy
- Wavelength setting repeatability
- Bandpass
- Baseline flatness
- Baseline stability
- Noise level



Example printout 1 of validation function with dot impact printer (for 1 check item)




Automatic calibration and self-diagnostic function


Upon initialization, memory check, wavelength drive mechanism check, lamp ignition check and 656.1nm check are automatically diagnosed and adjusted.

Also, the automatic wavelength calibration function is incorporated for automatic calibration by using the bright line of D₂ lamp as a reference wavelength.

Besides, each cumulative ignition time of the W1 and D₂ lamps can be displayed. This function is useful for instrument management.



Example printout 2 of validation function with dot impact printer (covering all check items)



Easy to Save Measured Data.

Measured data can be saved in a commercially available USB flash memory*.
The data is savable in the text file format.
A commercial spreadsheet program can processing the measured data.

* Not all of the commercial USB flash memory products are usable. Contact us separately.


Data saving with USB flash memory

Data can be transferred easily to PC by means of a USB flash memory.



Data loading

Because the measured data is saved in the text file format, it is sent directly to a commercial spreadsheet program. This function can be utilized for detailed data analysis, report generation, etc.



Personal computer-based instrument control and data processing.

The Model U-2900 becomes controllable from PC when the UV Solutions program (P/N 1J1-0210) is installed. Not only the U-2900 functions such as photometry and wavelength scan, but also data pasting to another application program, data transfer, report generation, etc. will be supported.


(The Model U-2910 is a special instrument for PC control.)

A PC compatible with Windows® 7 needs to be prepared separately.)




Simple to use

Buttons are laid out so as to follow the operating procedure. Hence, operations can be performed smoothly.




A rich variety of data processing functions

Spectrum can be expanded, contracted, smoothed, differentiated and integrated, and fundamental arithmetic calculations are applicable between spectra.




Powerful support for report generation

Measured data can be sent to Microsoft® Excel®. One click of the Report button suffices for data transmission.




A spectrum can be pasted to another application program.




Print preview

Before printing, the contents can be checked by the Print Preview command.




Measured data text convertible

Besides conversion into the ASCII text or J-CAMP format, spectra can be saved in Metafile format.



Windows® 7 and Microsoft® Excel® are trademarks of Microsoft Corporation, U.S.A. Microsoft® and Windows® are registered trademarks of U.S. Microsoft Corporation in USA and other countries.

Instrument application further expandable by UV Solutions program.



U-2910 is a special model for PC control.

Option package program (P/N 1J1-0211)

Capable of performing color calculation and color difference calculation.

Nucleic acid measurement program (P/N 1J1-0212)

Usable for checking extraction and refinement of nucleic acids such as DNA/RNA essential for genetic research. Data (photometric values at 230, 260 and 280 nm), and calculation results (260/280, nucleic acid concentration, protein concentration, molar concentration) can be displayed collectively on the screen.



GLP/GMP program (P/N 1J1-0213)

The burden of optical instrument performance check can be alleviated.

Item-dependent measurement mode and automatic measurement mode (when no sample is necessary) have also been prepared. In addition, a judgment tolerance range is settable so as to meet each analytical purpose.

Report generator program (P/N 1J1-0214)

Developed to enable freely customizing the report output format of measurement results. Each character size and position of report items and comment are settable and besides, mathematical calculations conventionally dependent on manual method can be executed by spreadsheet program functions.



Diverse optional accessories selectable for specific applications such as multiple samples and trace sample.

Auto sipper

(P/N 2J1-0100)

This computer-controlled sample sipper is provided with a sample recovery function and other versatile functions. In combination with an autosampler, this unit will further advance automation and labor saving in the preparatory stage. This sipper cannot control temperature.



Specifications

| | |
|-----------------------|----------------------------------|
| Minimum sample volume | 0.6 mL |
| Carryover | 1% or less |
| Cell capacity | About 50 μ L |
| Optical path length | 10 mm |
| Reference beam side | 10 mm rectangular cell mountable |

Auto sipper with temperature control

(P/N 2J1-0101)

Designed to have the same structure as the above auto sipper. Because the flow cell section can be maintained at a constant temperature, exact control is ensured.



Specifications (without temperature control on reference side)

| | |
|-----------------------|----------------------------------|
| Minimum sample volume | 0.6 mL |
| Carryover | 1% or less |
| Cell capacity | About 50 μ L |
| Optical path length | 10 mm |
| Settable temperature | 20 to 40°C |
| Reference beam side | 10 mm rectangular cell mountable |

AS-1010 Auto sampler

(P/N 2J1-0121/0122)

In combination with the auto sipper or in flow injection analysis, this unit is used for multiple-sample measurement. A suction needle can be moved in three directions X, Y and Z.

Specifications (sample tube not included)

| | |
|------------------|---|
| Sample tube size | Outer diameter 15 mm, height 105 mm (option required) |
| | Outer diameter 12 mm, height 105 mm |

6 cell positioner with temperature control

(P/N 2J1-0103/0104)

Six 10 mm cells can be mounted on the sample beam side, and they can be changed over automatically at certain intervals.



Specifications

| | |
|---|---|
| Difference in capacity due to cell changeover | Within $\pm 0.5\%$ (at 100%T) |
| Applicable cell | 10 mm cell (not included in this product) |
| Settable temperature | 20 to 40°C |

Electronic thermostatted cell holder

(P/N 131-0306/0307)

Standard-equipped with a magnet stirrer to maintain a constant sample temperature in each cell. Temperature can be indicated down to a minimum 0.1°C scale. Because of an electronic thermostatted type using forced air cooling, this cell holder is capable of quick heating and cooling. A thermostatic chamber is unnecessary. (Temperature control: S only)



Specifications

| | |
|------------------------------|--|
| Control temperature range | 10° to 60°C (settable in steps of 0.1°C) (normal temperature 25°C) |
| Temperature setting accuracy | Within $\pm 2^\circ\text{C}^*$ (difference between set temperature and sample temperature) |
| Temperature stability | Within $\pm 0.5^\circ\text{C}^*$ |

* Room temperature: 25°C, sample: distilled water

Electronic thermostatted cell holder

(Thermostatted water bath is needed separately)(P/N 131-0301/0302)

In the measurement of protein or nucleic acid fusion, a sample temperature can be changed continuously to determine variation in absorbance. Being of an electronic thermostatted type, this cell holder is capable of quick heating and cooling. Sample temperature can be increased and decreased isothermally. In addition, the set temperature can be maintained evenly inside a cell, because a stirrer is provided. (Temperature control: S and R)



Specifications

| | |
|--|--|
| Control temperature range | 0° to 100°C (settable in steps of 0.1°C) (normal temperature 25°C) |
| Temperature setting accuracy | Within $\pm 2^\circ\text{C}^*$ (difference between set temperature and sample temperature) |
| Temperature stability | Within $\pm 0.5^\circ\text{C}^*$ |
| Equipped with isothermal temperature regulation function | |

* Room temperature: 25°C, sample: distilled water, circulatory water temperature: 22°C, set temperature: -10 to 105°C
A circulating thermostatic chamber needs to be prepared separately.

Water circulated cell holder

(P/N 210-2111)

Water from a thermostatic chamber is circulated through this cell holder to maintain a sample cell at a constant temperature.



Specifications


(circulating thermostatic chamber and cell not included in this product)

| | |
|-----------------------------|---------------------------------|
| Operating temperature range | From normal temperature to 40°C |
| Temperature stability | Within $\pm 0.3^\circ\text{C}$ |

Micro flow cell unit

(P/N 210-2113)

Suitable for continuous measurement of a trace sample.



Specifications

| | |
|---------------------|---|
| Cell capacity | 70 μ L |
| Optical path length | 10 mm (quartz flow cell used) |
| Withstand pressure | Max. 0.1 MPa |
| Connecting tube | Fluoropolymer. Outer diameter 2 mm, inner diameter 1 mm |

Flow cell unit

(P/N 210-2173)

The inside of this cell is structured to minimize stagnation of liquid and adhesion of air bubbles.



Specifications

| | |
|---------------------|---|
| Cell capacity | 600 μ L |
| Optical path length | 5 mm (quartz flow cell used) |
| Withstand pressure | Max. 0.1 MPa |
| Connecting tube | Fluoropolymer. Outer diameter 4 mm, inner diameter 3 mm |
| Reference beam side | 5 mm rectangular cell (standard accessory) |

LC flow cell unit

(P/N 210-2131)

Exclusively used for a liquid chromatograph. In wavelength scan, the baseline is always kept flat.




Specifications

| | |
|---------------------|--|
| Cell capacity | 8 μ L |
| Optical path length | 8 mm (quartz flow cell used) |
| Baseline flatness | ± 0.001 Abs (200 to 350 nm) ± 0.002 Abs (190 to 850 nm) |

Micro cell holder

(P/N 122-0060)

Suitable for measurement of trace samples in medical and biochemical fields.



Specifications (cell required separately)

| | |
|--|--|
| Wavelength range | 220 to 950 nm |
| Cell mounting /dismounting repeatability | Within $\pm 0.3\%$ T |
| Baseline flatness | Within ± 0.005 Abs (50 mL micro-volume cell used) |

Tandem cell holder

(P/N 210-2115)

Up to three 10 mm cells can be mounted on each of the sample and reference beam sides. Sample temperature can be maintained at a constant level by circulating constant-temperature water through the cell holder section.



Specifications (circulating thermostatic chamber and cell not included in this product)

| | |
|---------------------------|-------------------------|
| Control temperature range | 15 to 40°C |
| Temperature stability | $\pm 0.3^\circ\text{C}$ |

Rectangular long path cell holder

(P/N 210-2107)




Conditions of applicable rectangular cell are given below.

| | |
|---------------------|----------------------------|
| Optical path length | 10, 20, 30, 40, 50, 100 mm |
| Outer width | 12.75 mm |

5 Position turretcell holder

(P/N 210-2110)

Five 10 mm rectangular cells can be mounted on the sample beam side and a micro cell mask (P/N 200-1537) can be inserted in each cell holder. (Cell and micro cell mask are not included in this product.) It is recommended to prepare a set of five cells.



Recommended cell sets

| | |
|----------|--|
| 124-0352 | 10 mm quartz cell set (5 cells in set) |
| 124-0378 | 10 mm glass cell set (5 cells in set) |

4 Position rectangular long path cell holder

(P/N 150-0940)

Four rectangular long path absorption cells can be mounted on the sample beam side and they can be changed over externally.



Specification

| | |
|-------------|---|
| Cell length | 100 mm, 50 mm to 10 mm cells applicable |
|-------------|---|

Glass filter holder

(P/N 210-2109)

Used for transmittance/absorbance measurement of such a solid sheet sample as glass filter.



Specifications

| | |
|------------------|---|
| Sample thickness | 0.5 to 5 mm |
| Sample size | 12 \times 25 mm to 55 \times 100 mm |

Cylindrical long path cell holder

(P/N 210-2108)

Used for measurement with a cylindrical cell.



Film holder

(P/N 210-2112)

Designed for measurement of film-like samples.



Specifications

| | |
|---------------|-----------------------------------|
| Film frame | Width 25 mm, height 30 to 50 mm |
| Beam aperture | Width 10 mm \times height 20 mm |

Polarizer holder

(P/N 210-2130)

Used for linearly polarizing the sample beam to check its polarization characteristic and for measuring optical rotary power with a sample placed between the polarizer and analyzer.



Specifications

| | |
|------------------|--|
| Wavelength range | 400 to 750 nm |
| Sample thickness | 0.5 to 5 mm |
| Sample size | Min. 12 \times 25 mm/max. 55 \times 100 mm |

Mask for Micro cell

To be inserted into a standard rectangular cell holder for measurement of a trace sample.

Specifications

| | |
|----------|------------------------------|
| 200-1537 | Mask for Micro cell (1.5 mm) |
| 200-1538 | Mask for Micro cell (1.2 mm) |

Cell

| | |
|----------|--------------------------------|
| 124-0357 | Micro quartz cell, 10 mm |
| 200-0551 | Black quartz micro cell, 10 mm |

The following cells are usable for the above micro cell holder (P/N 122-0060).



| Part No. | Part name | Capacity | Optical path length |
|----------|---------------|------------|---------------------|
| 130-0622 | 50 micro cell | 50 μ L | 10 mm |
| 130-0623 | 25 micro cell | 25 μ L | 5 mm |
| 130-0621 | 5 micro cell | 5 μ L | 0.5 mm |

Examined 10mm rectangular quartz cell

(P/N 210-1462)

Optical path lengths measured at 21 points on cell using a three-dimensional measuring instrument are indicated down to the fourth decimal place (mm).

Communication cable

(P/N 121-1521)

A computer and the spectrophotometer main unit need to be connected with this cable for computer control of measurement using the UV Solutions program.

U-2900/2910

Specifications

| | |
|---|--|
| Optical system | Double beam |
| Wavelength range | 190 to 1,100 nm |
| Spectral bandpass | 1.5 nm |
| Stray light | 0.05% or less (220 nm for NaI, 340 nm for NaNO ₂) |
| Wavelength accuracy | ±0.3 nm (at 656.1, 486.0 nm) |
| Wavelength setting repeatability | ±0.1 nm |
| Photometric range | –3 to 3 Abs 0 to 300%T |
| Photometric accuracy (certified according to NIST SRM 930) | ±0.002 Abs (0 to 0.5 Abs) |
| | ±0.004 Abs (0.5 to 1.0 Abs) |
| | ±0.008 Abs (1.0 to 2.0 Abs) |
| | ±0.3%T |
| Photometric repeatability (certified according to NIST SRM 930) | ±0.001 Abs (0 to 0.5 Abs) |
| | ±0.002 Abs (0.5 to 1.0 Abs) |
| | ±0.004 Abs (1.0 to 2.0 Abs) |
| | ±0.1%T |
| Wavelength scan speed | 10, 100, 200, 400, 800, 1,200, 2,400, 3,600 nm/min |
| Response | Fast, standard, slow |
| Baseline stability | 0.0003 Abs/h (at 500 nm, 2 hours after power-on) |
| Noise level | ±0.00015 Abs (at 500 nm) |
| Baseline flatness | ±0.0006 Abs (within 200 to 950 nm) |
| Light source | W1 and D ₂ lamps |
| Light source changeover | Auto (user selectable from 325 to 370 nm) |
| Detector | Silicon photodiode |
| Display | U-2900: color LCD with backlight (26.4 cm) |
| Printer I/F | U-2900: Centronics interface |
| Serial I/F | RS-232C (exclusive for UV Solutions program) |
| Size (main unit) | U-2900: 500 (W) × 605 (D) × 283 (H)mm (with LCD lowered) |
| | U-2910: 500 (W) × 605 (D) × 241 (H)mm (without PC and printer) |
| Weight (main unit) | U-2900: 31 kg, U-2910: 29 kg |
| Power supply | 100, 115, 220, 230 or 240 V, 50/60 Hz |
| Power consumption | 300 VA |

Software functions

- Measurement mode
 - Photometry
 - Wavelength scan
 - Time scan
 - Multiple-wavelength
- Ratio (260/280)
- Working curve type
 - Linear
 - Quadratic
 - Polygonal line
 - K factor input
- Calculation of correlation coefficient
- Concentration unit input
- Kinetic assay
- Spectrum and working curve printout
- Spectrum display
- Peak/valley detection
- Tracing
- Scale expansion/contraction
- Smoothing
- Differentiation
- Area calculation
- Fundamental arithmetic calculations between spectra
- Data saving
- Validation function
- Automatic wavelength calibration
- Lamp ignition time

* These are functions of Model U-2900 main unit. Equivalent functions are provided under PC control.

● Windows®, Microsoft®, Excel® are registered trademarks of Microsoft Corporation in the USA and other countries, and other company names and product names are registered trademarks, trademarks, or trade name of respective companies.



Science Ring

The Hitachi High-Tech Group aims to be a global leader in the "view", "measure" and "analyze" scientific and analysis fields, maintains points of contact with customers in a wide range of disciplines and actively works to provide advanced high added value solutions.

The logo mark is centered on the "S" from "Science", which represents the form created and connected through our cooperation as a good partner to customers and society that has its roots in our technologies, and which is expressed as organic spheres encircled by a ring.

It indicates our promise to society to create value through high-tech solutions that connect science and society.

CAUTION: For correct operation, follow the instruction manual when using the instrument.

Specifications in this catalog are subject to change with or without notice, as Hitachi High-Tech Science Corporation continues to develop the latest technologies and products for our customers.

NOTICE: The system is For Research Use Only, and is not intended for any animal or human therapeutic or diagnostic use. These data are an example of measurement; the individual values can not be guaranteed.

Hitachi High-Technologies Corporation

Tokyo, Japan

www.hitachi-hitec.com/global/science/

24-14, Nishi-shimbashi 1-chome, Minato-ku Tokyo 105-8717, Japan

For technical consultation before purchase, please contact:
contact@nst.hitachi-hitec.com

