

Tubes with Screw Cap



Polystyrene Tubes


Screw Cap, Conical / Round Bottom

- Further Closures p. 5 | 14
- Overview: See Technical Appendix for max. relative centrifugal force (RCF)
- Cell Culture Tubes p. 1 | 24



Cat.-No.	164 180	164 161	163 177
Bottom design	conical	conical	round
ø [mm] x height [mm]	16.8 x 100	16.8 x 100	17 x 100
Volume [ml]	12	12	12
Sterile	-	+	-
Screw cap	-	-	black
Bayonet cap ^{*)}	white	blue	-
Graduation	-	-	-
Writing area	-	-	-
Support skirt	+	+	-
Quantity per bag/case	1300	25/1000	900

^{*)} open by a 1/8 turn



Cat.-No.	188 161	188 171	186 161	186 171
Bottom design	conical	conical	round	round
ø [mm] x height [mm]	17 x 120	17 x 120	17 x 120	17 x 120
Volume [ml]	15	15	15	15
Sterile	+	+	+	+
Screw cap	blue	blue	blue	blue
Graduation	+	+	-	-
Writing area	+	+	-	-
Support skirt	-	-	-	-
Quantity per bag/case	50*/500	100/1000	50*/500	100/1000

^{*)} box

➤ New

➤ New

➤ New

➤ New

1 Cell / Tissue Culture
2 HTS- Microplates
3 Immunology / HLA
4 Microbiology / Bacteriology
5 Tubes / Multi-Purpose Beakers
6 Liquid Handling
7 Molecular Biology
8 Protein Crystallisation
9 Separation
10 Biochips / Microfluidics
11 Cryo-Technics
12 Lids / Sealers / Cap Mats
13 Reaction Tubes / Analyser Cups
14 Accessories




Polypropylene Tubes Screw Cap, Conical / Round Bottom

Overview: See Technical Appendix for max. relative centrifugal force (RCF)

- Very good mechanical, thermal and chemical stability
- Available in different packaging units and types
- All listed conical tubes are non-cytotoxic
- Light protection tubes for light-sensitive materials and reactions

DNase-free
RNase-free
human DNA-free
non-Pyrogenic




Cat.-No.	163 270	188 261	188 271	210 261	210 270	227 261
Bottom design	round	conical	conical	conical	conical	conical
ø [mm] x height [mm]	16.8 x 100	17 x 120	17 x 120	30 x 115	30 x 115	30 x 115
Volume [ml]	12	15	15	50	50	50
Sterile	-	+	+	+	+	+
Colour	natural	natural	natural	natural	natural	natural
Screw cap	white	blue	blue	blue	blue	blue
Graduation	-	blue	blue	blue	blue	blue
Writing area	-	white	white	white	white	white
Support skirt	-	-	-	+	+	-
Quantity per bag/case	900	50*/500	100/1000	25/450	25*/300	20/500

*) box

*) box

Light protection tubes



Cat.-No.	227 270	227 285	188 283	188 280	227 283	227 280
Bottom design	conical	conical	conical	conical	conical	conical
ø [mm] x height [mm]	30 x 115	30 x 115	17 x 120	17 x 120	30 x 115	30 x 115
Volume [ml]	50	50	15	15	50	50
Sterile	+	+	+	+	+	+
Colour	natural	natural	brown	brown	brown	brown
Screw cap	blue	white	blue	blue	blue	blue
Graduation	blue	blue	blue	blue	blue	blue
Writing area	white	white	white	white	white	white
Support skirt	-	-	-	-	-	-
Quantity per bag/case	25*/300	20/500	50*/500	100/1000	25*/300	20/500

*) box

➡ New

➡ New

➡ New

➡ New